

ZN644 rabbit pAb

Cat No.:ES19278

For research use only

Overview

Product Name	ZN644 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZN644 AA range: 674-724
Specificity	This antibody detects endogenous levels of human ZN644
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger protein 644 (Zinc finger motif enhancer-binding protein 2) (Zep-2)
Gene Name	ZN644 KIAA1221 ZEP2
Cellular localization	Nucleus.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84146
Human Swiss-Prot Number	Q9H582
Alternative Names	
Background	

