



SPAT5 rabbit pAb

Cat No.:ES19282

For research use only

Overview

Product Name	SPAT5 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SPAT5 AA range: 9-59
Specificity	This antibody detects endogenous levels of human SPAT5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Spermatogenesis-associated protein 5 (ATPase family protein 2 homolog) (Spermatogenesis-associated factor protein)
Gene Name	SPATA5 AFG2 SPAF
Cellular localization	Cytoplasm . Mitochondrion . Cytoplasm, cytoskeleton, spindle .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	166378
Human Swiss-Prot Number	Q8NB90
Alternative Names	
Background	

