



SE6L2 rabbit pAb

Cat No.:ES19288

For research use only

Overview

Product Name	SE6L2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SE6L2 AA range: 106-156
Specificity	This antibody detects endogenous levels of human SE6L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Seizure 6-like protein 2
Gene Name	SEZ6L2 PSK UNQ1903/PRO4349
Cellular localization	Cell membrane ; Single-pass type I membrane protein . Endoplasmic reticulum membrane ; Single-pass type I membrane protein . Detected on cell surface of lung-cancer.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	26470
Human Swiss-Prot Number	Q6UXD5
Alternative Names	
Background	

