

INT12 rabbit pAb

Cat No.:ES19290

For research use only

Overview

Product Name	INT12 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human INT12 AA range: 398-448
Specificity	This antibody detects endogenous levels of human INT12
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Integrator complex subunit 12 (Int12) (PHD finger protein 22)
Gene Name	INTS12 PHF22 SBB122
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	57117
Human Swiss-Prot Number	Q96CB8
Alternative Names	
Background	

