

# FA92B rabbit pAb

Cat No.:ES19292

For research use only

## Overview

Product Name	FA92B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FA92B AA range: 187-237
Specificity	This antibody detects endogenous levels of human FA92B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein FAM92B
Gene Name	FAM92B
Cellular localization	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Cytoplasm, cytoskeleton, cilium basal body . Extensive colocalization with CBY1 at mother centrioles. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	339145
Human Swiss-Prot Number	Q6ZTR7
Alternative Names	
Background	

