



LPP60 rabbit pAb

Cat No.:ES19293

For research use only

Overview

Product Name	LPP60 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LPP60 AA range: 112-162
Specificity	This antibody detects endogenous levels of human LPP60
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	60 kDa lysophospholipase (EC 3.1.1.5) [Includes: L-asparaginase (EC 3.5.1.1) (L-asparagine amidohydrolase); Platelet-activating factor acetylhydrolase (PAF acetylhydrolase) (EC 3.1.1.47)]
Gene Name	ASPG C14orf76
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	374569
Human Swiss-Prot Number	Q86U10
Alternative Names	
Background	

