

RN216 rabbit pAb

Cat No.: ES19294

For research use only

Overview

Product Name RN216 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse **Recommended dilutions** WB 1:1000-2000

Immunogen Synthesized peptide derived from human RN216 AA

range: 196-246

Specificity This antibody detects endogenous levels of human

RN216

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name E3 ubiquitin-protein ligase RNF216 (EC 6.3.2.-) (RING

finger protein 216) (Triad domain-containing protein 3) (Ubiquitin-conjugating enzyme 7-interacting protein 1) (Zinc finger protein inhibiting NF-

Gene Name RNF216 TRIAD3 UBCE7IP1 ZIN

Cellular localization Cytoplasm.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/ml

Observed band

Human Gene ID 54476 Human Swiss-Prot Number Q9NWF9

Alternative Names

+86-27-59760950

Background

