



ZWINT rabbit pAb

Cat No.:ES19301

For research use only

Overview

Product Name	ZWINT rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZWINT AA range: 165-215
Specificity	This antibody detects endogenous levels of human ZWINT
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ZW10 interactor (ZW10-interacting protein 1) (Zwint-1)
Gene Name	ZWINT
Cellular localization	Nucleus. Chromosome, centromere, kinetochore. Localizes to kinetochores from late prophase to anaphase.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	11130
Human Swiss-Prot Number	O95229
Alternative Names	
Background	

