



FNDC4 rabbit pAb

Cat No.:ES19303

For research use only

Overview

Product Name	FNDC4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FNDC4 AA range: 33-83
Specificity	This antibody detects endogenous levels of human FNDC4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Fibronectin type III domain-containing protein 4 (Fibronectin type III repeat-containing protein 1)
Gene Name	FNDC4 FRCP1 UNQ6389/PRO21134
Cellular localization	Membrane ; Single-pass type I membrane protein . Secreted . The N-terminus is probably cleaved to release a secreted extracellular portion of the protein. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	64838
Human Swiss-Prot Number	Q9H6D8
Alternative Names	
Background	

