



GBP6 rabbit pAb

Cat No.:ES19305

For research use only

Overview

Product Name	GBP6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human GBP6 AA range: 243-293
Specificity	This antibody detects endogenous levels of human GBP6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Guanylate-binding protein 6 (GTP-binding protein 6) (GBP-6) (Guanine nucleotide-binding protein 6)
Gene Name	GBP6
Cellular localization	Cytoplasmic vesicle .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	163351
Human Swiss-Prot Number	Q6ZN66
Alternative Names	
Background	

