

CAN14 rabbit pAb

Cat No.:ES19308

For research use only

Overview

Product Name	CAN14 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CAN14 AA range: 136-186
Specificity	This antibody detects endogenous levels of human CAN14
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Calpain-14 (EC 3.4.22.-) (Calcium-activated neutral proteinase 14) (CANP 14)
Gene Name	CAPN14
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	440854
Human Swiss-Prot Number	A8MX76
Alternative Names	
Background	

