



GNPTG rabbit pAb

Cat No.:ES19313

For research use only

Overview

Product Name	GNPTG rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human GNPTG AA range: 216-266
Specificity	This antibody detects endogenous levels of human GNPTG
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	N-acetylglucosamine-1-phosphotransferase subunit gamma (GlcNAc-1-phosphotransferase subunit gamma) (UDP-N-acetylglucosamine-1-phosphotransferase subunit gamma)
Gene Name	GNPTG C16orf27 GNPTAG CAB56184 LP2537
Cellular localization	Secreted . Golgi apparatus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84572
Human Swiss-Prot Number	Q9UJJ9
Alternative Names	
Background	

