

TDRD1 rabbit pAb

Cat No.:ES19330

For research use only

Overview

Product Name	TDRD1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TDRD1 AA range: 803-853
Specificity	This antibody detects endogenous levels of human TDRD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Tudor domain-containing protein 1 (Cancer/testis antigen 41.1) (CT41.1)
Gene Name	TDRD1
Cellular localization	Cytoplasm . Component of the meiotic nuage, also named P granule, a germ-cell-specific organelle required to repress transposon activity during meiosis. Also present in chromatoid body (By similarity). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	56165
Human Swiss-Prot Number	Q9BXT4
Alternative Names	
Background	

