

# PBX1 rabbit pAb

Cat No.:ES19338

For research use only

## Overview

Product Name	PBX1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PBX1 AA range: 162-212
Specificity	This antibody detects endogenous levels of human PBX1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Pre-B-cell leukemia transcription factor 1 (Homeobox protein PBX1) (Homeobox protein PRL)
Gene Name	PBX1 PRL
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	5087
Human Swiss-Prot Number	P40424
Alternative Names	
Background	

