

B3GN1 rabbit pAb

Cat No.:ES19345

For research use only

Overview

Product Name	B3GN1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human B3GN1 AA range: 13-63
Specificity	This antibody detects endogenous levels of human B3GN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	N-acetyllactosaminide beta-1,3-N-acetylglucosaminyltransferase (EC 2.4.1.149) (I-beta-1,3-N-acetylglucosaminyltransferase) (iGnT) (Poly-N-acetyllactosamine extension enzyme) (UDP-GlcNAc:betaGal beta-1
Gene Name	B3GNT1 B3GNT6
Cellular localization	Golgi apparatus membrane ; Single-pass type II membrane protein . Localizes near the trans-Golgi apparatus. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	11041
Human Swiss-Prot Number	O43505
Alternative Names	
Background	

