

PK1L2 rabbit pAb

Cat No.:ES19347

For research use only

Overview

Product Name	PK1L2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PK1L2 AA range: 1505-1555
Specificity	This antibody detects endogenous levels of human PK1L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Polycystic kidney disease protein 1-like 2 (PC1-like 2 protein) (Polycystin-1L2)
Gene Name	PKD1L2 KIAA1879 PC1L2
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	114780
Human Swiss-Prot Number	Q7Z442
Alternative Names	
Background	

