



D107A rabbit pAb

Cat No.:ES19348

For research use only

Overview

Product Name	D107A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human D107A AA range: 3-53
Specificity	This antibody detects endogenous levels of human D107A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Beta-defensin 107 (Beta-defensin 7) (BD-7) (DEFB-7) (Defensin, beta 107)
Gene Name	DEFB107A DEFB107 DEFB7; DEFB107B
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	245910
Human Swiss-Prot Number	Q8IZN7
Alternative Names	
Background	

