

IGSF1 rabbit pAb

Cat No.:ES19350

For research use only

Overview

Product Name	IGSF1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human IGSF1 AA range: 545-595
Specificity	This antibody detects endogenous levels of human IGSF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Immunoglobulin superfamily member 1 (IGSF1) (Immunoglobulin-like domain-containing protein 1) (Inhibin-binding protein) (InhBP) (Pituitary gland-specific factor 2) (p120)
Gene Name	IGSF1 IGDC1 KIAA0364 PGSF2
Cellular localization	[Isoform 1]: Membrane ; Multi-pass membrane protein .; [Isoform 2]: Membrane ; Multi-pass membrane protein .; [Isoform 3]: Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3547
Human Swiss-Prot Number	Q8N6C5
Alternative Names	
Background	

