

# SAC31 rabbit pAb

Cat No.:ES19354

For research use only

## Overview

Product Name	SAC31 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SAC31 AA range: 31-81
Specificity	This antibody detects endogenous levels of human SAC31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	SAC3 domain-containing protein 1 (SAC3 homology domain-containing protein 1)
Gene Name	SAC3D1 SHD1
Cellular localization	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cytoplasm, cytoskeleton, spindle . Localizes on centrosomes in interphase cells and at spindles in mitosis. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	A6NKF1
Alternative Names	
Background	

