



NTNG2 rabbit pAb

Cat No.:ES19360

For research use only

Overview

Product Name	NTNG2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NTNG2 AA range: 384-434
Specificity	This antibody detects endogenous levels of human NTNG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Netrin-G2 (Laminin-2)
Gene Name	NTNG2 KIAA1857 LMNT2 UNQ9381/PRO34206
Cellular localization	Cell membrane ; Lipid-anchor, GPI-anchor ; Extracellular side .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84628
Human Swiss-Prot Number	Q96CW9
Alternative Names	
Background	

