



# MFS11 rabbit pAb

Cat No.:ES19365

For research use only

## Overview

Product Name	MFS11 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MFS11 AA range: 323-373
Specificity	This antibody detects endogenous levels of human MFS11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	UNC93-like protein MFSD11 (Major facilitator superfamily domain-containing protein 11) (Protein ET)
Gene Name	MFSD11 ET
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79157
Human Swiss-Prot Number	O43934
Alternative Names	
Background	

