

KV124 rabbit pAb

Cat No.:ES19367

For research use only

Overview

Product Name	KV124 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KV124 AA range: 5-55
Specificity	This antibody detects endogenous levels of human KV124
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Ig kappa chain V-I region Daudi
Gene Name	
Cellular localization	Secreted . Cell membrane .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	P04432
Alternative Names	
Background	

