



SYCC rabbit pAb

Cat No.:ES19369

For research use only

Overview

Product Name	SYCC rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SYCC AA range: 532-582
Specificity	This antibody detects endogenous levels of human SYCC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Cysteine--tRNA ligase, cytoplasmic (EC 6.1.1.16) (CysteinyI-tRNA synthetase) (CysRS)
Gene Name	CARS
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	833
Human Swiss-Prot Number	P49589
Alternative Names	
Background	

