



PRAF1 rabbit pAb

Cat No.:ES19374

For research use only

Overview

Product Name	PRAF1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PRAF1 AA range: 130-180
Specificity	This antibody detects endogenous levels of human PRAF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Prenylated Rab acceptor protein 1 (PRA1 family protein 1)
Gene Name	RABAC1 PRA1 PRAF1
Cellular localization	Cell membrane ; Multi-pass membrane protein . Cytoplasm . Golgi apparatus . Cytoplasmic vesicle, secretory vesicle, synaptic vesicle . According to some authors, it is an integral membrane protein, while others showed that it is cytoplasmic and membrane-associated to Golgi and synaptic vesicles. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10567
Human Swiss-Prot Number	Q9UI14
Alternative Names	
Background	

