



FSCN2 rabbit pAb

Cat No.:ES19375

For research use only

Overview

Product Name	FSCN2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FSCN2 AA range: 356-406
Specificity	This antibody detects endogenous levels of human FSCN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Fascin-2 (Retinal fascin)
Gene Name	FSCN2
Cellular localization	Cytoplasm, cytoskeleton . Cell projection, stereocilium .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	25794
Human Swiss-Prot Number	O14926
Alternative Names	
Background	

