

# CASA1 rabbit pAb

Cat No.:ES19379

For research use only

## Overview

Product Name	CASA1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CASA1 AA range: 45-95
Specificity	This antibody detects endogenous levels of human CASA1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Alpha-S1-casein [Cleaved into: Casoxin-D]
Gene Name	CSN1S1 CASA CSN1
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	1446
Human Swiss-Prot Number	P47710
Alternative Names	
Background	

