

RASN rabbit pAb

Cat No.:ES19384

For research use only

Overview

Product Name	RASN rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RASN AA range: 7-57
Specificity	This antibody detects endogenous levels of human RASN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	GTPase NRas (Transforming protein N-Ras)
Gene Name	NRAS HRAS1
Cellular localization	Cell membrane ; Lipid-anchor ; Cytoplasmic side . Golgi apparatus membrane ; Lipid-anchor . Shuttles between the plasma membrane and the Golgi apparatus. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	4893
Human Swiss-Prot Number	P01111
Alternative Names	
Background	

