



TRAT1 rabbit pAb

Cat No.:ES19386

For research use only

Overview

Product Name	TRAT1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TRAT1 AA range: 56-106
Specificity	This antibody detects endogenous levels of human TRAT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	T-cell receptor-associated transmembrane adapter 1 (T-cell receptor-interacting molecule) (TRIM) (pp29/30)
Gene Name	TRAT1 TCRIM HSPC062
Cellular localization	Cell membrane ; Single-pass type III membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	50852
Human Swiss-Prot Number	Q6PIZ9
Alternative Names	
Background	

