

CYTSB rabbit pAb

Cat No.:ES19387

For research use only

Overview

Product Name	CYTSB rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CYTSB AA range: 403-453
Specificity	This antibody detects endogenous levels of human CYTSB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Cytospin-B (Nuclear structure protein 5) (NSP5) (Sperm antigen HCMOGT-1) (Sperm antigen with calponin homology and coiled-coil domains 1)
Gene Name	SPECC1 CYTSB NSP5
Cellular localization	Nucleus .; [Isoform 3]: Membrane ; Lipid-anchor .; [Isoform 4]: Membrane ; Lipid-anchor .; [Isoform 5]: Membrane ; Lipid-anchor .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	92521
Human Swiss-Prot Number	Q5M775
Alternative Names	
Background	

