



CN080 rabbit pAb

Cat No.:ES19390

For research use only

Overview

Product Name	CN080 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CN080 AA range: 31-81
Specificity	This antibody detects endogenous levels of human CN080
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Uncharacterized protein C14orf80
Gene Name	C14orf80
Cellular localization	Cell projection, cilium . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	283643
Human Swiss-Prot Number	Q86SX3
Alternative Names	
Background	

