



TEX36 rabbit pAb

Cat No.:ES19391

For research use only

Overview

Product Name	TEX36 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TEX36 AA range: 102-152
Specificity	This antibody detects endogenous levels of human TEX36
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Testis-expressed sequence 36 protein
Gene Name	TEX36 C10orf122
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	387718
Human Swiss-Prot Number	Q5VZQ5
Alternative Names	
Background	

