

# CM035 rabbit pAb

Cat No.:ES19392

For research use only

## Overview

Product Name	CM035 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CM035 AA range: 2-52
Specificity	This antibody detects endogenous levels of human CM035
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Putative uncharacterized protein C13orf35
Gene Name	C13orf35
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	400165
Human Swiss-Prot Number	Q6ZP68
Alternative Names	
Background	

