

EVI2A rabbit pAb

Cat No.:ES19393

For research use only

Overview

Product Name	EVI2A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EVI2A AA range: 135-185
Specificity	This antibody detects endogenous levels of human EVI2A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Protein EVI2A (Ecotropic viral integration site 2A protein homolog) (EVI-2A)
Gene Name	EVI2A EVDA EVI2
Cellular localization	Membrane; Single-pass type I membrane protein.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	2123
Human Swiss-Prot Number	P22794
Alternative Names	
Background	

