

SP9 rabbit pAb

Cat No.:ES19398

For research use only

Overview

Product Name	SP9 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SP9 AA range: 187-237
Specificity	This antibody detects endogenous levels of human SP9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	Transcription factor Sp9
Gene Name	SP9
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	100131390
Human Swiss-Prot Number	POCG40
Alternative Names	
Background	

