

TRH rabbit pAb

Cat No.:ES19399

For research use only

Overview

Product Name	TRH rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TRH AA range: 133-183
Specificity	This antibody detects endogenous levels of human TRH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Prothyroliberin [Cleaved into: Thyroliberin (Protirelin) (TSH-releasing factor) (Thyrotropin-releasing factor) (TRF) (Thyrotropin-releasing hormone) (TRH)]
Gene Name	TRH
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7200
Human Swiss-Prot Number	P20396
Alternative Names	
Background	

