

UD2B7 rabbit pAb

Cat No.:ES19403

For research use only

Overview

Product Name	UD2B7 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human UD2B7 AA range: 467-517
Specificity	This antibody detects endogenous levels of human UD2B7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	UDP-glucuronosyltransferase 2B7 (UDPGT 2B7) (EC 2.4.1.17) (3,4-catechol estrogen-specific UDPGT) (UDP-glucuronosyltransferase 2B9) (UDPGT 2B9) (UDPGTh-2)
Gene Name	UGT2B7 UGTB2B9
Cellular localization	Endoplasmic reticulum membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7364
Human Swiss-Prot Number	P16662
Alternative Names	
Background	

