

# KRA31 rabbit pAb

Cat No.:ES19407

For research use only

## Overview

Product Name	KRA31 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KRA31 AA range: 11-61
Specificity	This antibody detects endogenous levels of human KRA31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Keratin-associated protein 3-1 (High sulfur keratin-associated protein 3.1) (Keratin-associated protein 3.1)
Gene Name	KRTAP3-1 KAP3.1 KRTAP3.1
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	83896
Human Swiss-Prot Number	Q9BYR8
Alternative Names	
Background	

