



XKR2 rabbit pAb

Cat No.:ES19409

For research use only

Overview

Product Name	XKR2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human XKR2 AA range: 376-426
Specificity	This antibody detects endogenous levels of human XKR2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	XK-related protein 2 (Membrane protein XPLAC) (X Kell blood group-related, X-linked)
Gene Name	XKRX XKR2 XPLAC XRG2
Cellular localization	Cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	402415
Human Swiss-Prot Number	Q6PP77
Alternative Names	
Background	

