

# COX7S rabbit pAb

Cat No.:ES19410

For research use only

## Overview

Product Name	COX7S rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human COX7S AA range: 13-63
Specificity	This antibody detects endogenous levels of human COX7S
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Putative cytochrome c oxidase subunit 7A3, mitochondrial (Cytochrome c oxidase subunit VIIa 3)
Gene Name	COX7A2P2 COX7A3 COX7AL2 COX7AP2
Cellular localization	Mitochondrion inner membrane .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	O60397
Alternative Names	
Background	

