



TIM44 rabbit pAb

Cat No.:ES19412

For research use only

Overview

Product Name	TIM44 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TIM44 AA range: 260-310
Specificity	This antibody detects endogenous levels of human TIM44
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Mitochondrial import inner membrane translocase subunit TIM44
Gene Name	TIMM44 MIMT44 TIM44
Cellular localization	Mitochondrion inner membrane ; Peripheral membrane protein ; Matrix side . Mitochondrion matrix .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10469
Human Swiss-Prot Number	O43615
Alternative Names	
Background	

