



TR16L rabbit pAb

Cat No.:ES19413

For research use only

Overview

Product Name	TR16L rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TR16L AA range: 91-141
Specificity	This antibody detects endogenous levels of human TR16L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Tripartite motif-containing protein 16-like protein (Tripartite motif-containing protein 70)
Gene Name	TRIM16L TRIM70
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	147166
Human Swiss-Prot Number	Q309B1
Alternative Names	
Background	

