

CF089 rabbit pAb

Cat No.:ES19414

For research use only

Overview

Product Name	CF089 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CF089 AA range: 193-243
Specificity	This antibody detects endogenous levels of human CF089
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Bombesin receptor-activated protein C6orf89 (Amfion)
Gene Name	C6orf89 BRAP UNQ177/PRO203
Cellular localization	[Isoform 1]: Golgi apparatus membrane ; Single-pass type II membrane protein . Midbody .; [Isoform 2]: Golgi apparatus membrane ; Single-pass type II membrane protein . Midbody .; [Isoform 3]: Cytoplasm . Nucleus, nucleolus . Retained in nucleolar organiser regions (NORs) in mitotic cells (PubMed:23460338). .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	221477
Human Swiss-Prot Number	Q6UWU4
Alternative Names	
Background	

