



CLU rabbit pAb

Cat No.:ES19416

For research use only

Overview

| | |
|---------------------------------|--|
| Product Name | CLU rabbit pAb |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human;Mouse |
| Recommended dilutions | WB 1:1000-2000 |
| Immunogen | Synthesized peptide derived from human CLU AA range: 435-485 |
| Specificity | This antibody detects endogenous levels of human CLU |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C . Avoid repeated freeze-thaw cycles. |
| Protein Name | Clustered mitochondria protein homolog |
| Gene Name | KIAA0664 |
| Cellular localization | Cytoplasm . Cytoplasmic granule . A fraction colocalizes with tyrosinated tubulin and can be detected close to mitochondria. . |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | 23277 |
| Human Swiss-Prot Number | O75153 |
| Alternative Names | |
| Background | |

