



DUS23 rabbit pAb

Cat No.:ES19424

For research use only

Overview

Product Name	DUS23 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DUS23 AA range: 67-117
Specificity	This antibody detects endogenous levels of human DUS23
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Dual specificity protein phosphatase 23 (EC 3.1.3.16) (EC 3.1.3.48) (Low molecular mass dual specificity phosphatase 3) (LDP-3) (VH1-like phosphatase Z)
Gene Name	DUSP23 LDP3 VHZ
Cellular localization	Cytoplasm, cytosol. Nucleus. Mainly cytosolic. Also nuclear.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	54935
Human Swiss-Prot Number	Q9BVJ7
Alternative Names	
Background	

