



OD3L2 rabbit pAb

Cat No.:ES19425

For research use only

Overview

Product Name	OD3L2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human OD3L2 AA range: 52-102
Specificity	This antibody detects endogenous levels of human OD3L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Outer dense fiber protein 3-like protein 2
Gene Name	ODF3L2 C19orf19
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	284451
Human Swiss-Prot Number	Q3SX64
Alternative Names	
Background	

