



DSA2D rabbit pAb

Cat No.:ES19427

For research use only

Overview

Product Name	DSA2D rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DSA2D AA range: 38-88
Specificity	This antibody detects endogenous levels of human DSA2D
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Putative dispanin subfamily A member 2d (DSPA2d)
Gene Name	
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	C9JQL5
Alternative Names	
Background	

