



NFU1 rabbit pAb

Cat No.:ES19428

For research use only

Overview

Product Name	NFU1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NFU1 AA range: 203-253
Specificity	This antibody detects endogenous levels of human NFU1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	NFU1 iron-sulfur cluster scaffold homolog, mitochondrial (HIRA-interacting protein 5)
Gene Name	NFU1 HIRIP5 CGI-33
Cellular localization	Mitochondrion. Cytoplasm, cytosol.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	27247
Human Swiss-Prot Number	Q9UMS0
Alternative Names	
Background	

