



LNX1 rabbit pAb

Cat No.:ES19432

For research use only

Overview

Product Name	LNX1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LNX1 AA range: 527-577
Specificity	This antibody detects endogenous levels of human LNX1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	E3 ubiquitin-protein ligase LNX (EC 6.3.2.-) (Ligand of Numb-protein X 1) (Numb-binding protein 1) (PDZ domain-containing RING finger protein 2)
Gene Name	LNX1 LNX PDZRN2 UNQ574/PRO1136
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84708
Human Swiss-Prot Number	Q8TBB1
Alternative Names	
Background	

