



ALKB8 rabbit pAb

Cat No.:ES19433

For research use only

Overview

Product Name	ALKB8 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ALKB8 AA range: 27-77
Specificity	This antibody detects endogenous levels of human ALKB8
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Alkylated DNA repair protein alkB homolog 8 (EC 1.14.11.-) (Probable alpha-ketoglutarate-dependent dioxygenase ABH8) (S-adenosyl-L-methionine-dependent tRNA methyltransferase ABH8) (tRNA (carboxymethyl) ALKBH8 ABH8
Gene Name	ALKBH8 ABH8
Cellular localization	Cytoplasm . Nucleus . Predominantly cytoplasmic. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	91801
Human Swiss-Prot Number	Q96BT7
Alternative Names	
Background	

