



PAR3L rabbit pAb

Cat No.:ES19434

For research use only

Overview

Product Name	PAR3L rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PAR3L AA range: 587-637
Specificity	This antibody detects endogenous levels of human PAR3L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Partitioning defective 3 homolog B (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 19 protein) (PAR3-beta) (Partitioning defective 3-like protein) (PAR3-L protein)
Gene Name	PARD3B ALS2CR19 PAR3B PAR3L
Cellular localization	Endomembrane system. Cell junction. Cell junction, tight junction. Partially localized along the cell-cell contact region. Colocalizes with TJP1 to epithelial tight junctions.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	117583
Human Swiss-Prot Number	Q8TEW8
Alternative Names	
Background	

